


- ▶ Tracheostomy tube holder - 4 variants: “K”, “H”, “K PED”, “H PED”
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Consists of a light foam core, 3 cm wide, which is surrounded by a soft non-woven fabric
- ▶ Thin adherent strips, which are threaded into the eyelets of the tracheostomy tube’s neck flange
- ▶ White color
- ▶ “K” variant: with adherent fastener
- ▶ “H” variant: with hook fastener
- ▶ “K PED” variant: for pediatric use, with adherent fastener
- ▶ “H PED” variant: for pediatric use, with hook fastener
- ▶ Individually packed



## Order information:

<b>SUPRAFIX® 11 H PED, 15 - 26 cm</b>	<b>REF 32308</b>
<b>SUPRAFIX® 11 K PED, 18 - 25 cm</b>	<b>REF 32509</b>
<b>SUPRAFIX® 13 H PED, 17 - 28 cm</b>	<b>REF 32310</b>
<b>SUPRAFIX® 13 K PED, 20 - 27 cm</b>	<b>REF 32510</b>
<b>SUPRAFIX® 15 H PED, 19 - 30 cm</b>	<b>REF 32311</b>
<b>SUPRAFIX® 15 K PED, 22 - 29 cm</b>	<b>REF 32512</b>
<b>SUPRAFIX® 18 H, 22 - 41 cm</b>	<b>REF 32301</b>
<b>SUPRAFIX® 18 K, 25 - 40 cm</b>	<b>REF 32502</b>
<b>SUPRAFIX® 25 H, 35 - 42 cm</b>	<b>REF 32300</b>
<b>SUPRAFIX® 25 K, 32 - 48 cm</b>	<b>REF 32500</b>
<b>SUPRAFIX® 30 H, 34 - 54 cm</b>	<b>REF 32305</b>
<b>SUPRAFIX® 30 K, 37 - 53 cm</b>	<b>REF 32501</b>
<b>SUPRAFIX® 36 H, 40 - 58 cm</b>	<b>REF 32303</b>
<b>SUPRAFIX® 36 K, 43 - 57 cm</b>	<b>REF 32503</b>



Single patient use



Consult instructions for use



CE marking